


UNIVERSITÀ
DEGLI STUDI
DI PADOVA




UNIVERSITÀ
DEGLI STUDI
DI PADOVA

Acquisizione delle sequenze in risonanza magnetica cardiaca Parte 1


Anna Baritussio, MD, PhD, FEACVI




A “one stop shop” technique



BIVENTRICULAR
VOLUMES AND FUNCTION




MYOCARDIAL OEDEMA




T1-WEIGHTING





LATE ENHANCEMENT



INDUCIBLE ISCHEMIA



T1/T2 MAPPING



VALVE ASSESSMENT



RESEARCH

Open Access

Standardized cardiovascular magnetic resonance imaging (CMR) protocols: 2020 update



Schulz-Menger et al. *Journal of Cardiovascular Magnetic Resonance* (2020) 22:19
<https://doi.org/10.1186/s12968-020-00610-6>

Christopher M. Kramer^{1*} , Jörg Barkhausen², Chiara Bucciarelli-Ducci³, Scott D. Flamm⁴, Raymond J. Kim⁵ and Eike Nagel⁶

POSITION STATEMENT


Open Access

Standardized image interpretation and post-processing in cardiovascular magnetic resonance - 2020 update



Society for Cardiovascular Magnetic Resonance (SCMR): Board of Trustees Task Force on Standardized Post-Processing

Jeanette Schulz-Menger^{1*}, David A. Bluemke², Jens Bremerich³, Scott D. Flamm⁴, Mark A. Fogel⁵, Matthias G. Friedrich⁶, Raymond J. Kim⁷, Florian von Knobelsdorff-Brenkenhoff⁸ , Christopher M. Kramer⁹, Dudley J. Pennell¹⁰, Sven Plein¹¹ and Eike Nagel¹²




SAGITTAL




CORONAL

Free breathing



AXIAL

Bright blood



Black blood



Rodrigues *et al. Journal of Cardiovascular Magnetic Resonance* (2016) 18:26
DOI 10.1186/s12968-016-0246-1

Journal of Cardiovascular
Magnetic Resonance

REVIEW

Open Access




Extra-cardiac findings in cardiovascular magnetic resonance: what the imaging cardiologist needs to know

Jonathan C. L. Rodrigues^{1,2,3}, Stephen M. Lyen², William Loughborough², Antonio Matteo Amadu^{1,4}, Anna Baritussio^{1,5}, Amardeep Ghosh Dastidar^{1,5}, Nathan E. Manghat² and Chiara Bucciarelli-Ducci^{1,5*}


Meta-analysis - 12 studies, 7062 patients

- Prevalence extra-cardiac findings 35% (major findings in 12%)
- Change in management in 1%


May contribute to diagnosis
of cardiac pathology




Right thyroid lobe enlargement with local mass effect on the trachea




Right thyroid lobe enlargement with local mass effect on the trachea




Incidental florid mediastinal adenopathy




Right thyroid lobe enlargement with local mass effect on the trachea




Incidental florid mediastinal adenopathy



Extra-medullary haematopoiesis



**LUPV draining into
left brachiocephalic vein**





LV SV (ml/mq)

RV SV (ml/mq)

67


149



Absent pericardium



Lung Interposition
between Aorta and Pulmonary Artery




Levoposition
Cardiac Apex Pointing Posteriorly



Image property	Spoiled Gradient Echo	bSSFP
Blood to myocardium contrast	Variable – relies on through plane flow	Good – T_2/T_1 ratio
Flow sensitivity	High – can visualise flow jets	Low – uses balanced gradients
Image quality	Low SNR – no shimming required	High SNR – good shimming required


Cine images are acquired during a breath-hold.

Breath-hold on expiration provides more consistent positioning but inspiratory breath-hold may be more comfortable and easier to sustain for some patients.




Transaxial

Bisecting
Mitral Valve
and Apex




Pseudo 2-chamber




Pseudo 2-chamber


Bisecting
Mitral Valve
and Apex








Pseudo 4-chamber



True short axis view














AVOIDING LVOT




AVOIDING RV










Biventricular volumes and function assessment
free of any geometric assumption






Cine Images




Best identification of tricuspid valve plane




Cine Images



4C



3C




4CH




Sagittal RVOT

★ 3-point planning



AXIAL



RV in-/outflow





Image property	Spoiled Gradient Echo	bSSFP
Blood to myocardium contrast	Variable – relies on through plane flow	Good – T_2/T_1 ratio
Flow sensitivity	High – can visualise flow jets	Low – uses balanced gradients
Image quality	Low SNR – no shimming required	High SNR – good shimming required


Cine images are acquired during a breath-hold.

Breath-hold on expiration provides more consistent positioning but inspiratory breath-hold may be more comfortable and easier to sustain for some patients.




Gradient echo cine sequences







MITRAL REGURGITATION




AORTIC STENOSIS




LVOT OBSTRUCTION



MITRAL REGURGITATION



AORTIC STENOSIS




LVOT OBSTRUCTION

N.B. a perfect 2D image plane
of a 3D structure is impossible



- High spatial resolution
- High contrast blood-cusps

Table 1. CMR Pulse Sequences for Valvular Heart Disease

Pulse Sequence	Indication
<u>SSFP cine</u>	Valve anatomy and motion Ventricular volumes and function
Gradient echo cine	Valve anatomy and motion Turbulent flow
Phase contrast	Velocity Forward and regurgitant volumes
Turbo spin echo (T1, T2, with or without fat saturation)	Valve mass characterization
Segmented inversion recovery gradient echo	Valve mass characterization



Mitral Stenosis



Phase-contrast sequences

“Bright blood” GRE



Acquired on a plane perpendicular to flow direction

Level of grey depending on flow velocity and direction

- still protons → grey


- protons moving perpendicular to the plane → black or white

VENEC – systemic flow 150 cm/s, pulmonary flow 100 cm/s




Valve
stenosis/regur-
gitation,
 $Q_p:Q_s$,
sistemic-
pulmonary
shunt


Phase-contrast sequences




Sagittal LVOT




Coronal LVOT




Ao flow



Sagittal RVOT



Coronal RVOT




MPA flow


VENC settings:

- Optimal within 25% of the true peak velocity
 - Too low: flow aliasing
 - Too high: underestimating velocity
- Correct direction of flow (Right-Left, Foot-Head)
- Image plane distal from valve leaflet tips


Phase-contrast sequences






Sagittal LVOT




Coronal LVOT




Ao flow




Sagittal RVOT




Coronal RVOT




MPA flow




Phase-contrast sequences




Sagittal LVOT




Coronal LVOT




Ao flow





Sagittal RVOT




Coronal RVOT



MPA flow








Aortic Assessment



1. Aortic annulus
2. Sinuses of Valsalva
3. Sinotubular junction
4. Mid ascending aorta (level of the pulmonary arteries)
5. Proximal aortic arch (at the origin of the brachiocephalic trunk)
6. Mid aortic arch (btw. left common carotid and subclavian arteries)
7. Proximal descending thoracic aorta
8. Mid descending aorta (level of the pulmonary arteries)
9. At diaphragm
10. At the celiac axis origin
11. Immediately proximal to aortic bifurcation.


MR Angiography










MR Angiography



Partial anomalous pulmonary venous return -
right superior lobar vein draining into the SVC





UNIVERSITÀ
DEGLI STUDI
DI PADOVA